



Zettachring 6, 70567 Stuttgart, Germany
 Tel: +49-711-13257-0 Fax: +49-711-13257-89
 E-Mail: info.sep.eu@tosoh.com
 Web: www.tosohbioscience.com

3604 Horizon Drive, Suite 100, King of Prussia, PA 19406, USA
 Tel: +1 800-366-4875 Fax: +1 610-272-3028
 E-Mail: info.sep.am@tosoh.com
 Web: www.tosohbioscience.com

OPERATING CONDITIONS and SPECIFICATIONS

TSK-GEL® SP-2SW Products

Part Numbers: 07165 4.6 mm ID x 25.0 cm L SP-2SW 5 µm

Functional Group: -C₃H₆-SO₃Na

Small Ion Capacity > 0.3 meq/g

This sheet contains the recommended operating conditions and the specifications for TSK-GEL SP-2SW columns. Installation instructions and column care information is described in a separate Instruction Manual.

A. OPERATING CONDITIONS

1. Shipping Solvent: Methanol
 Replace the shipping solvent with water before changing to the mobile phase of choice

2. Max. Flow Rate: 1.2 ml/min 5.0 cm L column
 1.0 ml/min 25.0 cm L column

NOTE When a buffer with high viscosity is used, the maximum flow rate may have to be reduced so as not to exceed the maximum pressure drop. When changing solvents, use a flow rate equal to 50% of the maximum flow rate.

3. Standard Flow Rate: 0.6 - 0.8 ml/min 25.0 cm L

4. Max. Pressure: 150 kg/cm² = 2200 psi 25.0 cm L

5. pH Range 2.0 - 7.5

6. Ionic Strength Typical 0.5 M - 0.1 M starting buffer and 0.2 M and 0.5 M final buffer.

7. Organic Conc. 0 - 100% for aqueous soluble organic solvents. Make gradual solvent changes using a shallow gradient at low flow rate.

8. Temperature: 10 - 45°C. Reduce flow rate when operating below 10°C.

9. Cleaning Solvents: (1) Buffer with 0.5 - 1.0 M salt, or
 (2) pH 2 - 3 buffer, or
 (3) Buffer containing methanol or acetonitrile, or, if nothing else was successful,
 (4) Urea or non-ionic surfactant in buffer

NOTE Choose a cleaning solvent based on sample properties, e.g. use
 (1) to remove strongly adsorbed proteins,
 (2) to reduce the negative charge on the biopolymer,
 (3) to decrease hydrophobic adsorption, and
 (4) to increase protein solubility.

10. Storage: Store the column in methanol at ambient temperature when it will not be used the next day. For overnight storage flush the column with the low salt concentration mobile phase at 0.2 ml/min. Prevent air from entering the column!

11. Column Protection The use of guard columns is recommended to prolong the life of the analytical column. Guard column life depends greatly on sample cleanliness. As a general rule, guard columns should be replaced when the peaks become excessively wide, or when the peaks show splitting.

B. SPECIFICATIONS

The performance of TSK-GEL SP-2SW columns is tested under the conditions described in the Data Sheet. All columns have passed the following quality control specifications:

Number of Theoretical Plates (N): ≥ 5,000 25.0 cm L
 Asymmetry Factor (AF): 0.8 - 1.6 25.0 cm L

Note our technical hotline tel +49 711 13257-57
 and e-mail, techsupport.sep@tosoh.com