




TOSOH



# Sensitive Quantification of AAVs and their Impurities



Analysis of AAVs

## Your Challenge

- ▶ You need a method that detects AAVs at low titers.
- ▶ You need to separate AAVs from impurities.

## Our Solution

LenS<sub>3</sub><sup>™</sup> MALS detector & TSKgel<sup>®</sup> GMPW<sub>XL</sub> column

- ▶ Sensitive detector and reliable separations

What was done?

- ▶ AAV5 and AAV8 samples were separated by size exclusion chromatography and detected by MALS

What was the result?

- ▶ Limit of detection for AAV5 is  $8 \times 10^9$  vc & host-cell impurities are efficiently separated from AAVs

**Determine down to  $8 \times 10^9$  AAV capsids and separate impurities of AAV samples by combining the efficient separation by a mixed bed SEC column and the highly sensitive MALS detector LenS<sub>3</sub>.**

## Your Benefit

**Reduce sample volume and detect lowest AAV titers**


TOSOH BIOSCIENCE



# SEC-MALS Sensitively Quantifies Adeno-Associated Viruses and Detects Impurities

Cell and gene therapy in medicine relies on efficient and specific delivery of therapeutic genetic material to target cells. These delivery devices are called vectors, and many have been developed from viruses<sup>1</sup>. The parvovirus family of viruses called adeno-associated viruses (AAVs), have emerged as the vector of choice in therapeutic gene delivery due to their optimized tissue and cell specificity and their ability to edit the genome precisely. AAVs are non-enveloped viruses with an icosahedral structure belonging to the parvovirus family and containing single-stranded DNA (ssDNA) (~4.7 kb in length) as their genome (*Figure 1*). To date, 13 AAV serotypes and >100 variants have been identified with variable tissue tropism<sup>2-4</sup>.

**Figure 1.** Graphic representation of a wild-type AAV and its genome.



The viral ssDNA cassette contains replicating genes (Rep) for virus multiplication, capsid genes (Cap), and poly(A) tail (pA). Cap encodes three structural viral proteins VP1, VP2, and VP3, needed for the capsid assembly<sup>3,4</sup>.

The figure was created using BioRender.com from the original source: <https://www.dynotx.com/introduction-to-aav-as-a-gene-therapy-vector-part-1/>

Development of clinically desirable AAV capsids with optimal genome design requires rapid, accurate, and robust analytical methods to assess AAV purity, capsid titer, DNA content, and structure-activity relationships<sup>5</sup>. Size exclusion chromatography (SEC), when coupled with multi-angle light scattering (MALS), offers a powerful analytical method for AAV characterization. In contrast to the traditional method of AAV analysis using the UV 260/280 nm absorbance ratio, the advantage of MALS detection is the high sensitivity for AAVs due to their considerable capsid mass (~3.7 MDa), inherently providing a strong light scattering response for analytical characterization. Here, we describe the use of SEC-MALS to quantify AAV capsid titers and to determine sample impurities.

## Experimental Conditions

### Samples

Purified AAV5 ( $2.0 \times 10^{13}$  viral capsids (vc)/mL) was purchased from Virovek (Hayward, CA). AAV8-containing human embryonic kidney (HEK) cell culture supernatant (titer  $1.5 \times 10^{12}$  vc/mL) was from ExcellGene SA (Monthey, Switzerland). Bovine serum albumin (BSA) standard was purchased from Thermo Scientific (Rockford, IL).

### SEC-MALS Analysis

AAV samples were analyzed on a TSKgel GMPW<sub>XL</sub> SEC column coupled with the LenS<sub>3</sub> MALS detector. Capsid titers during analyses were as follows: AAV5 at  $2.0 \times 10^{11}$  vc/mL was used to compare light scattering signals (right, low, and high angles) to UV absorbance at 280 nm. For LOQ/LOD determination, AAV5 dilutions from  $1.0 \times 10^{12}$  vc/mL to  $1.6 \times 10^9$  vc/mL were used. AAV8-containing supernatant was tested at  $1.5 \times 10^{12}$  vc/mL and  $1.5 \times 10^{11}$  vc/mL titers. The SEC-MALS system was calibrated with BSA prior to sample analysis to determine the detectors' offsets, band-broadening parameters, and response factors. All data were processed and analyzed using the SECview software.

Instrument: Thermo Fisher Scientific UltiMate™  
3000 U/HPLC  
Column: TSKgel GMPW<sub>XL</sub> (7.8 mm ID × 30 cm,  
13 μm particles, pores 125-1000 Å)  
Detectors: (1) Ultimate 3000 variable wavelength  
detector (VWD) at 280 nm  
(2) Tosoh LenS<sub>3</sub> MALS detector  
Mobile phase: 50 mmol/L HEPES, 150 mmol/L NaCl, pH 7.1  
Flow rate: 0.5 mL/min  
Temperature: 25 °C  
Injection vol.: 40-50 μL


## Results and Discussion

### SEC-MALS Analysis

For the 22-25 nm sized AAV capsids, a TSKgel GMPW<sub>XL</sub> analytical column was chosen for separation due to its mixed pore size distribution up to 1000 Å. This SEC column enables excellent resolution for the analysis of AAVs, and their separation from accompanying impurities (cell culture media components, aggregated or degraded capsids, etc.).

Figure 2 shows the right-angle light scattering (RALS) signals for AAV5 and BSA. Characterization of AAV5 by SEC-MALS using the TSKgel GMPW<sub>XL</sub> column revealed a major AAV5 monomer peak, with a known capsid size of ~3.7 MDa for AAVs6, eluting at 14-15 min retention time. Some unassembled capsid proteins and other components were detected at ~20 min retention time. BSA (66 kDa), which was used for system calibration, eluted significantly later than AAV5 due to its much smaller size (insert in Figure 2).

Figure 2. SEC-MALS detection (RALS) of AAV5 and BSA.



**MALS vs. UV Absorbance at 280 nm**

Due to the high molecular weight (MW) of AAVs, the MW-dependent response of light scattering provided much higher sensitivity than the regular UV absorbance at 280 nm, as illustrated in Figure 3. When analyzing AAV5 at  $2 \times 10^{11}$  vc/mL (40  $\mu$ L injection,  $8 \times 10^9$  vc total loaded), no UV signal was detected at 280 nm, whereas strong MALS signals (signal/noise ratio of 10 to 60) were obtained at three different light scattering angles: RALS (90°), LALS (10°), HALS (170°).

Figure 3. Comparison of RALS, LALS, and HALS signals to UV @ 280 nm for AAV5.





Figure 4. Limit of detection by RALS signal.



**Linearity and Sensitivity**

The linearity and sensitivity of AAV detection by MALS were characterized by determining the limit of detection/limit of quantitation (LOD/LOQ) using serial dilutions from the starting AAV5 concentration of  $1.0 \times 10^{12}$  vc/mL and loading 40  $\mu$ L injections onto the TSKgel GMPW<sub>XL</sub> column (Figure 4).


The left panel shows 40  $\mu\text{L}$  injections of serially diluted samples from  $1.0 \times 10^{12}$  vc/mL ( $4.0 \times 10^{10}$  vc loaded) to  $1.6 \times 10^9$  vc/mL ( $6.4 \times 10^7$  vc loaded). The right panel shows the linearity of RALS signal from the same dilutions. The insert in the right panel magnifies the plot with the lowest AAV5 concentrations. As low as  $7.0 \times 10^9$  vc/mL AAV5 samples ( $2.8 \times 10^8$  vc total loaded) were detectable by SEC-MALS using the LenS<sub>3</sub> MALS detector's signal (RALS). The LOQ was then calculated from the LOD to be  $2.3 \times 10^{10}$  vc/mL ( $9.2 \times 10^8$  vc total loaded). Additionally, RALS detection on the LenS<sub>3</sub> provided a strong linear response across a large concentration range, allowing the determination of capsid titers.

### AAV8-containing Cell Culture Supernatant

An additional challenge of AAV analysis is to achieve efficient separation of AAVs from the remaining impurities to allow the quantification of AAV capsids. To test the separation performance, we analyzed AAV8 directly from the HEK cell culture media without prior AAV8 purification using TSKgel GMPW<sub>XL</sub> column.

**Figure 5** shows the chromatographic separation of AAV8 capsids with the elution peak at  $\sim 14.5$  min retention time whereas the HEK cell culture impurities mainly eluted at 16-20 min retention time, with additional tailing at higher loading concentration. The injection of two different AAV8 capsid loads (50  $\mu\text{L}$  of  $1.5 \times 10^{12}$  vc/mL and 50  $\mu\text{L}$  of  $1.5 \times 10^{11}$  vc/mL) illustrates the concentration-dependent response of the RALS signal, allowing AAV separation and quantification from more complex sample matrices. AAV8 at both concentrations shows a clear peak whereas cell culture impurities saturated the RALS signal (upper part not shown).

**Figure 5.** SEC-MALS analysis of AAV8-containing cell culture supernatant.



## Conclusion

A rapid and reproducible analytical methodology is crucial for the timely development of high-quality tools for cell and gene therapy applications. This application note describes a powerful analytical SEC-MALS technique for the sensitive detection of AAVs using the LenS<sub>3</sub> MALS detector and the TSKgel GMPW<sub>XL</sub> column. The limit of detection (LOD) for the method was determined to be as low as  $7.0 \times 10^9$  vc/mL ( $2.8 \times 10^8$  vc loaded in 40  $\mu\text{L}$ ). The high sensitivity of SEC-MALS detection allows the injection of low-concentration AAV samples which is a highly sought attribute in the early development stages with often limited availability of material. In addition, AAV separation directly from HEK cell culture supernatant was possible enabling direct MALS detection without prior AAV purification. In summary, the SEC-MALS methodology described here provides a sensitive, powerful, and robust tool for the detection of AAVs during product development, viral vector production, and throughout quality control. Analytical characterization (virus size, molecular weight, empty/full ratio) based on SEC-MALS data will be fully described in a separate application note.

## References

- Zhao, Z, Anselmo, A.C, Mitragotri, S. Viral vector-based gene therapies in the clinic. *Bioeng. Transl. Med.* 7, e10258. DOI:10.1002/btm2.10258 (2022)
- Bijlani, S, Pang, K.M, Sivanandam, V, Singh, A, Chatterjee, S. The role of recombinant AAV in precise genome editing. *Front. Genome Ed.* 3, 799722. DOI=10.3389/fgeed.2021.799722 (2022)
- Pupo, A, Fernández, A, Low, S.H, François, A, Suárez-Amarán, L, Samulski, R.J. AAV vectors: The Rubik's cube of human gene therapy. *Mol. Ther.* 30: 3515-3541. DOI: 10.1016/j.ymthe.2022.09.015 (2022)
- Kang, L, Jin, S, Wang, J, Lv, Z, Xin, C, Tan, C, Zhao, M, Wang, L, Liu, J. AAV vectors applied to the treatment of CNS disorders: Clinical status and challenges. *J. Cont. Release* 355: 458-473. DOI: 10.1016/j.jconrel.2023.01.067 (2023)
- Gimpel, A.L, Katsikis, G, Sha, S, Maloney, A.J, Hong, M.S, Nguyen, T.N.T, Wolfrum, J, Springs, S.L, Sinsky, A.J, Manalis, S.R, Barone, P.W, Braatz, R.D. Analytical methods for process and product characterization of recombinant adeno-associated virus-based gene therapies. *Mol. Therapy: Methods Clin. Dev.* 20: 740-754. DOI: 10.1016/j.omtm.2021.02.010 (2021)
- Serrano, M.A.C, Furman, R, Chen, G, Tao, L. Mass spectrometry in gene therapy: Challenges and opportunities for AAV analysis. *Drug Discov. Today* 28:103442. DOI: 10.1016/j.drudis.2022.103442 (2023)

Tosoh Bioscience and TSKgel are registered trademarks of Tosoh Corporation. LenS is a trademark of Tosoh Bioscience LLC. UltiMate is a trademark of Thermo Fisher Scientific.

Featured Products		
0008025	TSKgel GMPW <sub>XL</sub>	7.8 mm ID $\times$ 30 cm, 13 $\mu\text{m}$ , 12.5-100 nm pore size
0040000	LenS <sub>3</sub> MALS detector	