



mAb PURIFICATION PLATFORM DEVELOPMENT USING SkillPak™ PRE-PACKED COLUMNS

INTRODUCTION

Hybridoma technology is an efficient method for production of monoclonal antibodies (mAbs). After eliciting an immune response in mice, the antibody-producing B-cells are harvested before fusion of IgG-expressing B cells with myeloma cells to allow propagation in cell culture.

This note details a case study where a proprietary mouse mAb was efficiently purified using a two-step purification process from hybridoma cell line supernatant (“feedstock”). The challenge was to obtain a few milligrams of mAb for cell-based functional assays to evaluate the purification process in spite of very low mAb titer (~0.03 g/L) in the feedstock. The two-step purification platform was established using SkillPak™ 1 mL and 5 mL columns from Tosoh Bioscience.

MATERIALS AND METHODS

Pre-packed columns

The SkillPak 1 mL and 5 mL columns (specifications listed in Table 1 and 2 respectively) are designed for fast method development or resin screening. They are pre-packed with TOYOPEARL®, TSKgel® or Ca⁺⁺Pure-HA® process chromatography media for separation of biomolecules such as monoclonal antibodies, proteins, and oligonucleotides. These columns are designed to be operated with commonly used low or medium pressure liquid chromatography systems.

SPECIFICATIONS OF SkillPak 1 mL COLUMN

| | |
|----------------------------|---|
| Column dimension | 7 mm ID × 2.5 cm bed height |
| Volume | 1 mL |
| Maximum flow rate | 4 mL/min (600 cm/hr) |
| Maximum operating pressure | 0.3 MPa |
| Connections | Standard fittings (10-32 for 1/16 inch capillary) |
| Shipping buffer | 20% ethanol for TOYOPEARL and TSKgel, 20 mmol/L phosphate with 20% ethanol for Ca ⁺⁺ Pure-HA |

➤ Table 1

SPECIFICATIONS OF SkillPak 5 mL COLUMN

| | |
|--------------------------------------|--|
| Column dimension | 8 mm ID × 10 cm bed height |
| Volume | 5 mL |
| Standard flow rate | 1.3 mL/min (150 cm/hr) |
| Maximum flow rate | 5 mL/min (600 cm/hr) for TOYOPEARL M and C grade resins; 2.5 mL/min (300 cm/hr) for TOYOPEARL S-grade, TSKgel, Ca ⁺⁺ Pure-HA and TOYOPEARL F grade resins |
| Maximum operating pressure | 0.3 MPa for TOYOPEARL resins, ≤ 0.4 MPa for TSKgel resins and Ca ⁺⁺ Pure-HA |
| Connections | Standard fittings (10-32 for 1/16 inch capillary) |
| Shipping buffer | 20% ethanol for TOYOPEARL and TSKgel, 20 mmol/L phosphate with 20% ethanol for Ca ⁺⁺ Pure-HA |
| Asymmetry factor (As) specifications | 0.8-1.4 for TOYOPEARL and TSKgel, 0.8-2.6 for Ca ⁺⁺ Pure-HA |

➤ Table 2

The study used a SkillPak 5 mL column pre-packed with TOYOPEARL AF-rProtein A HC-650F, a polymethacrylate, high dynamic binding capacity protein A affinity media, in the capture step and a SkillPak 1 mL column pre-packed with TOYOPEARL Sulfate-650F, a high salt-tolerant cation exchange (CEX) media, in the polishing step.

STEP 1 (CAPTURE). COLUMN CHARACTERISTICS

| | |
|----------------------------|---|
| Media | TOYOPEARL AF-rProtein A HC-650F |
| Column | SkillPak 5 mL TOYOPEARL AF-rProtein A HC-650F |
| Bed size | 8 mm ID × 10 cm |
| Particle size | 45 µm |
| Pore diameter | 100 nm |
| DBC* (5 min) | 70 g/L |
| DBC* (2 min) | 50 g/L |
| Caustic stability | >200 CIP cycles (0.1 mol/L NaOH) |
| Maximum operating pressure | 0.3 MPa |

➤ Table 3

STEP 2 (POLISH). COLUMN CHARACTERISTICS

| | |
|----------------------------|--------------------------------------|
| Media | TOYOPEARL Sulfate-650F |
| Column | SkillPak 1 mL TOYOPEARL Sulfate-650F |
| Bed size | 7 mm ID x 2.5 cm |
| Particle size | 45 µm |
| Pore diameter | 100 nm |
| DBC | >120 g/L of IgG |
| Caustic stability | >200 CIP cycles (0.1 mol/L NaOH) |
| Maximum operating pressure | 0.3 MPa |

Table 4

MONOCLONAL ANTIBODY

Proprietary mouse mAb

PURIFICATION PROTOCOL

Determination of hybridoma mAb elution pH on protein A

Column: SkillPak 5 mL TOYOPEARL AF-rProtein A HC-650F

Equilibration buffer: 0.1 mol/L Na₂HPO₄/NaH₂PO₄, 0.15 mol/L NaCl, pH 7.3

Mobile phase A (gradient): 0.1 mol/L acetate (NaOH), pH 4.5

Mobile phase B (gradient): 0.1 mol/L acetic acid, pH 2.9

Elution gradient: linear from pH 4.5 to 2.9 over 10 CV

Flow (load): 150 cm/hr (1.25 mL/min), 4 min residence time

Flow (wash/gradient): 240 cm/hr (2.0 mL/min) (ÄKTA™ avant 25 instrument)

Sample: 45 mL hybridoma cell culture supernatant (buffer-adjusted)

Capture method for purification of hybridoma mAb

Column: SkillPak 5 mL TOYOPEARL AF-rProtein A HC-650F

Equilibration buffer: 0.1 mol/L Na₂HPO₄/NaH₂PO₄, 0.15 mol/L NaCl, pH 7.3

Post-loading 1st wash: equilibration buffer (5 CV)

Post-loading 2nd wash: 0.1 mol/L acetate (NaOH), pH 5.0 (5 CV)

Elution: 0.1 mol/L acetate (NaOH), pH 3.5 (5 CV)

Column strip: 0.1 mol/L acetic acid, pH 2.9 (3 CV)

Column cleaning: 0.2 mol/L NaOH (3 CV), 15 min hold

Flow (load): 150 cm/hr (1.25 mL/min), 4 min residence time

Flow (wash/elution): 240 cm/hr (2.0 mL/min) (ÄKTA avant 25)

Temperature: ambient (room temperature)

Sample: 340 mL hybridoma cell culture supernatant (buffer-adjusted)

Polishing method to remove mAb aggregates and impurities

Column: SkillPak 1 mL TOYOPEARL Sulfate-650F

Equilibration buffer: 0.1 mol/L acetate (NaOH), pH 5.0

Post-loading 1st wash: equilibration buffer (5 CV)

Post-loading 2nd wash: 0.1 mol/L acetate (NaOH), 0.1 mol/L NaCl, pH 5.0 (10 CV)

Elution: 0.1 mol/L acetate (NaOH), 0.35 mol/L NaCl, pH 5.0 (10 CV)

Column strip: 0.1 mol/L acetate (NaOH), 1.0 mol/L NaCl, pH 5.0 (7 CV)

Flow (all steps): 156 cm/hr (1.0 mL/min) 4 min residence time (ÄKTA avant 25 instrument)

Sample: 21 mL (0.41 mg total mAb) diluted protein A eluate, (9.8mg/4.8mL)/5 = 0,41 mg/mL

Temperature: ambient (room temperature)

DETERMINATION OF mAb ELUTION PH ON TOYOPEARL AF-rProtein A HC-650F MEDIA




Figure 1

mAb CAPTURE ON TOYOPEARL AF-rProtein A HC-650F MEDIA




Figure 2

RESULTS AND DISCUSSIONS

Finding the appropriate buffer and pH for mAb elution reduces the risk of increasing aggregation in the mAb sample. To identify optimized conditions for mAb binding and elution, a buffer-adjusted mAb-containing hybridoma cell line supernatant was loaded onto a SkillPak 5 mL column pre-packed with TOYOPEARL AF-rProtein A HC-650F media. Figure 1 shows that a sharp mAb peak with elution max at pH 4.0 was obtained using a linear pH gradient. To maximize recovery, pH 3.5 was selected for step elution.

For prep scale purification, 340 mL of hybridoma cell line supernatant (titer ~0.03 g/L) was loaded on the SkillPak 5 mL column. After a short wash at pH 5.0, a sharp and efficient elution peak was obtained at the start of the pH 3.5 elution (Figure 2).

Total mAb recovery was 9.8 mg in the elution peak (4.8 mL). Eluate was prepared for a polishing step on a SkillPak 1 mL column packed with TOYOPEARL Sulfate-650F by adding four eluate volumes of 0.2 mol/L acetate, pH 5.0.

Aggregate and other impurities in the protein A eluate were removed by using the SkillPak 1 mL TOYOPEARL Sulfate-650F column. **Figure 3** shows an efficient elution of mAb at a 0.35 mol/L NaCl step in the equilibration buffer. Impurities were removed by the polishing step as can be seen before and after the elution of the mAb peak.

POLISHING CHROMATOGRAPHY ON TOYOPEARL SULFATE-650F




Figure 3

The mAb eluate from the cation exchange step was subjected to size exclusion chromatography (SEC) to confirm purity and the monomeric state of the collected mAb. The elution peak at ~7.5 minutes indicated a largely monomeric mAb with >98% purity (Figure 4).

SEC ANALYSIS OF mAb PURITY




Figure 4

CONCLUSION

The two-step process presented here using SkillPak pre-packed columns was effective for purification of this mAb from hybridoma cell line supernatant and the process can feasibly be scaled up to a pilot plant and eventually to manufacturing scale.

TOYOPEARL, TSKgel and Tosoh Bioscience are registered trademarks of Tosoh Corporation.
 Ca**Pure-HA is a registered trademark of Tosoh Bioscience LLC.
 SkillPak is a trademark of Tosoh Bioscience LLC.
 ÄKTA is a trademark of GE Healthcare Bio-Sciences.