



**RESIN INFORMATION SHEET**

**Product Name** TOYOPEARL® AF-rProtein L-650F  
(Functionalised resin for affinity chromatography)

**Part Numbers**

0023486	TOYOPEARL AF-rProtein L-650F, 10 mL
0023487	TOYOPEARL AF-rProtein L-650F, 25 mL
0023488	TOYOPEARL AF-rProtein L-650F, 100 mL
0023489	TOYOPEARL AF-rProtein L-650F, 1 L
0023490	TOYOPEARL AF-rProtein L-650F, 5 L

**Product Description** TOYOPEARL chromatographic resins are based on a rigid methacrylic polymer, resulting in high mechanical and chemical stability. Resins are available as non-functionalized “HW” series resins for size exclusion separations, and derivatized with surface chemistries for alternative modes of chromatography such as ion exchange, hydrophobic interaction or affinity separations.

TOYOPEARL<sup>®</sup> AF-rProtein L-650F chromatographic resin is a functionalized support resin for affinity chromatography. The resin is prepared by multipoint attachment of a recombinant Protein L ligand onto TOYOPEARL HW-65F. TOYOPEARL AF-rProtein L-650F is useful for the chromatographic purification of immunoglobulins and their fragments at lab and process scale. The rProtein L ligand is not animal derived.

<b>Operating Conditions</b>	Packing pressure	Typically 0.3 MPa
	Shipping solvent	20 % (v/v) ethanol
	Shipping formulation	72 % (v/v) slurry in shipping solvent (*)
	Pressure limiting factor	Depend on column hardware (typically 0.7 MPa)
	Operating linear flowrate	Typically 10 - 600 cm/hour
	Long-term storage conditions	20 % (v/v) ethanol at 2 - 8 °C
	Cleaning-in-place/Sanitization	0.03 - 0.1 mol/L NaOH

<b>Specifications</b>	Particle size distribution (min. 80 % within range)	30 - 60 µm
	Protein adsorption capacity (of human antibody)	Min. 64 g/L - gel
	Bacterial count	Max. 100 CFU/mL
	Endotoxin concentration	Max. 10.0 EU/mL
	Eluable matter	Max. 0.2 % (for dry gel)
	Foreign substance (colored particle)	Max. 6

<b>Additional Information</b>	Appearance	White resin slurry which settles upon standing
	Mean pore diameter (base resin)	100 nm (*)

(\*) The value is for reference only, not guaranteed.

Lot-specific data are included in the Certificate of Analysis (COA) shipped with the product. For detailed test procedures please refer to the appropriate Regulatory Support File.