

**RESIN INFORMATION SHEET**

Product Name TOYOPEARL® AF-Tresyl-650M
(Activated base resin for affinity chromatography)

Part Numbers

0014471	TOYOPEARL AF-Tresyl-650M, 5 g
0014472	TOYOPEARL AF-Tresyl-650M, 100 g
0014905	TOYOPEARL AF-Tresyl-650M, 200 g

Product Description TOYOPEARL chromatographic resins are based on a rigid methacrylic polymer, resulting in high mechanical and chemical stability. Resins are available as non-functionalized “HW” series resins for size exclusion separations, and derivatized with surface chemistries for alternative modes of chromatography such as ion exchange, hydrophobic interaction or affinity separations.

TOYOPEARL AF-Tresyl-650M is an activated support resin for affinity chromatography. The product is prepared by introducing tresyl groups onto TOYOPEARL HW-65. Tresyl-activated material is used to immobilize ligands with free amino or thiol groups, quickly and with high efficiency.

Operating Conditions	Packing pressure	Typically 0.3 MPa
	Shipping preservative	Dry powder
	Pressure limiting factor	Depend on column hardware (typically 0.7 MPa)
	Operating linear flowrate	Typically 10 - 600 cm/hour
	Long-term storage conditions	Dry at 2 - 8 °C
Specifications	Particle size distribution (min. 80 % within range)	40 - 90 µm
	Protein binding capacity	Not less than 60 mg/g - gel
	Foreign substance (colored particle)	Unobserved
Additional Information	Appearance	White powder
	Mean pore diameter (base resin)	100 nm (*)

(*) The value is for reference only, not guaranteed.

Lot-specific data are included in the Certificate of Analysis (COA) shipped with the product.
For detailed test procedures please refer to the appropriate Regulatory Support File.